



Fetoly™ — A Diagnoly technology


# Your Trusted Real-Time **AI Ultrasound Assistant**

Smart. Simple. Secure.

Even today, up to **50%** of fetal abnormalities remain undetected during pregnancy **and we are changing that.**



Constantly detects +100  
fetal anatomical  
structures



Live measurements  
supporting pathology  
detection



Auto capture and  
annotation of images



CE 0197

FETOLY V2026.02.25.1389

# The Fetoly application helps clinicians detect fetal abnormalities and verify exam quality.

## 1 DOWNLOAD

Get started by downloading the Fetoly app on your tablet.

## 2 PLUG

Link your device to your ultrasound system via HDMI

## 3 RUN

Recieve real-time AI guidance without changing your scanning routine



No cloud or internet needed



Integration with your PACS system



Works with major ultrasound systems



REQUEST A FREE TRIAL

✉ [contact@diagnoly.com](mailto:contact@diagnoly.com)

🌐 [www.diagnoly.com](http://www.diagnoly.com)